


Microwave Heating: Theory and Practice




Based on Lecture Prepared by Dr. Lois B. Jassie

Characteristics of Microwave Energy

- ◆ Microwaves are electromagnetic radiation between the far IR and radio waves
- ◆ Microwaves are nominally between 1 mm and 100 cm in length (*e.g.*, 2450 MHz wave is 12.25 cm)
- ◆ Microwave energy is non-ionizing, low photon electromagnetic radiation at the powers used
- ◆ Microwave radiation causes molecular (particle) and ionic motion and dipole rotation
- ◆ Microwave energy does not cause a change in molecular structure


CEM

Electromagnetic Spectrum



CEM

A Microwave



E = electric field
 H = magnetic field
 λ = wavelength (12.2 cm for 2450 MHz)
 c = speed of light (300,000 km/s, 3.00×10^8 m/s)

CEM

Microwave Energy *versus* Other Electromagnetic Energy

Radiation Type	Typical Frequency (MHz)	Quantum Energy (ev)	Chemical Bond Type	Chemical Bond Energy (ev)
Gamma Rays	3.0×10^{14}	1.24×10^6	H-OH	5.2
X-Rays	3.0×10^{13}	1.24×10^5	H-CH ₃	4.5
Ultraviolet	1.0×10^9	4.1	H-NHCH ₃	4.0
Visible Light	6.0×10^8	2.5	H ₃ C-CH ₃	3.8
Infrared Light	3.0×10^6	0.012	PhCH ₂ -COOH	2.4
Microwaves	2450	1.013×10^{-5}	H-O-H	0.21
Radio	1	4×10^{-9}		

CEM

Material Interaction With Microwaves

Conductor **Reflective**
 Metals reflect microwave energy and do not heat.

Insulator **Transparent**
 Numerous materials are transparent to microwave energy and will not heat but can be good insulators.

Dielectric **Absorptive**
 These materials absorb microwave energy and are heated.

CEM

Liquids (mineral acids, solvents) heat rapidly when exposed to microwave energy. Absorption of microwave energy occurs by two mechanism:



Dipole Rotation Ionic Conduction

CEM


Sources of (Internal) Heat (via Energy Transformations)

- ◆ Molecular Rotation
- ◆ Conformational Changes
- ◆ 3-Dimensional Distortion
- ◆ Ion Flow Enhancements
- ◆ Liquid Structure Dissipation

CEM



Schematic of Sample Heating by Microwaves



Reprinted with permission from Neas, E.; Collins, M. in "Introduction to Microwave Sample Preparation: Theory and Practice," copyright 1988 by The American Chemical Society.

CEM

Has this ever happened to you?

You try to make you coffee in the morning by heating the water in the microwave and when you add the coffee it boils all over the Counter top?

Explanation {Superheated}:
Under microwave heating water boils at 105°C until nucleated to boil at 100°C.



© H. M. Skip Kingston

Superheated Temperatures of Solvents Irradiated with Microwave Energy

Solvent	Boiling Point, °C **	Superheated Temperature, °C	Temperature Difference, °C
Water	100	105	5
1-butanol	117	138	21
2-butanol	98	127	29
tert-butanol	83	112	29
methanol	65	84	19
2-propanol	82	108	26
1-pentanol	136	157	21
2-pentanol	119	135	16
tert-pentanol	120	115	13
1-heptanol	176	208	32
ethylene glycol	196	216	20
acetone	56	89	33
ethyl acetate	77	102	25
chloroform	61	89	28
diethyl ether	35	60	25
tetrahydrofuran (THF)	67	103	36
acetonitrile	82	120	38
cyclohexane	155	186	31
methyl ethyl ketone (MEK)	80	110	30

Reference: Majetch, G.; Neas, E.; Hoopes, T. *Journal of Chemical Education*, 1994
 * In 1 liter flask; ** B.P. at 760 mm. Handbook of Chemistry and Physics, CRC Press, Inc.



Factors That Influence Solution Heating

Physical Properties of the Solution

Viscosity
 Temperature
 Polarity
 Heat Capacity
 Dielectric

Ion Characteristics (conduction only)

Concentration
 Charge
 Size
 Mobility

Wavelength




Factors that Affect Materials Heated with Microwaves

- ◆ Angle of incidence of radiation
- ◆ Frequency
- ◆ Dielectric Constant
- ◆ Impedance
- ◆ Loss Mechanisms
- ◆ Mass and Molecular Size
- ◆ Magnetic Properties



Variation of Penetration with Frequency for Water at 25°C



Dipole Moment

CCl ₄	0
CO	0.10
HCl	1.08
H ₂ S	1.10
C ₃ H ₈ O	1.66
HF	1.82
H ₂ O	1.85
CH ₃ Cl	1.87
HNO ₃	2.17
HCN	2.93
CsF	7.87

Handbook of Chemistry and Physics, 60th Edition, CRC Press, 1980.
 Table E-66*, Gas phase molecule



Dissipation Factor Tangent delta

Ratio of the sample's dielectric loss (loss factor), E''
 to its dielectric constant, E'

$$\text{Tangent } \delta = E''/E'$$

E'' Loss factor is a measure of the sample material's ability to convert electromagnetic energy to thermal energy (heat).

E' Dielectric constant is the sample material's ability to store the microwave energy.



Effect of Increasing NaCl Concentration on the Dissipation Factor

Molal Concentration	Tangent δ ($\times 10^{-4}$)
0.0	1570
0.1	2400
0.3	4350
0.5	6250

Reprinted with permission from Neas, E.; *Introduction to Microwave Sample Preparation: Theory and Practice*, copyright 1988, American Chemical Society, Washington, DC; Ch. 2. Measurements made at 3000 MHz, and 25 °C.

CEM

Effect of Temperature on the Dissipation Factor of Water

Temperature, °C	Tangent δ (10^{-4}) ^a
1.5	3100
5.0	2750
15.0	2050
25.0	1570
35.0	1270
45.0	1060
55.0	890
65.0	765
75.0	660
85.0	547
95.0	470

^a measurement at 3000 MHz and 25° C

CEM

Selected Physical and Dielectric Constants of Organic Solvents

Solvent	BP	VP	E'	Dipole Moment	tan $\delta \times 10^{-4}$
acetone	40	436	8.93	1.14	
methanol	56	184	20.7	2.69	6400
tetrahydrofuran	65	125	32.7	2.87	
hexane	69	120	1.88	<0.1	
ethyl acetate	77	73	6.02	1.88	
ethanol	78	-	24.3	1.69	2500
acetonitrile	82	89	37.5	3.44	
2-propanol	82	32	19.9	1.66	6700
1-propanol	97	14	20.3	3.09	~2400*
iso-octane	99	49	1.94	0	
water	100	760	78.3	1.87	1570
MIBK	116	20	13.11	-	
DMF	153	2.7	36.71	3.86	
DMSO	189	0.6	46.68	3.1	
ethylene glycol	198	-	41.0	2.3	10,000

note: data from Burdick & Jackson Laboratories Solvent Handbook. BP at 760 torr, VP (torr) at 25°C; e' at 20 °C; dipole moment at 25 °C. * at 10 °C (E.Peterson, IMPI symposium, 1989). Tan δ values from von Hippel, MIT Press, 1954.

CEM

Thermal and Microwave Characteristics of Laboratory Container Materials

Material	Melting Point (°C)	Maximum Service Temperature (°C)	Tangent δ ($\times 10^{-4}$)
Water			1570
Sodium chloride (0.1 molal)			2400
Polysulfone	<190	160	760
Phenol/formaldehyde	dec	120-190	519
Bakelite (asbestos filled)	dec	200-218	438
Nylon 6/6	253	102	128
Glass (Corning 0800)	>1000	---	126
Glass (Borosilicate)	>1080	---	12-75
Ceramic (depends on type)	---	---	6-50
Polypropylene	168-171	100-105	57
Polymethylmethacrylate	115	76-88	57
Porcelain (4462)	---	---	11
Polystyrene	242	82-91	3.3
Polyethylene	120-135	71-93	3.1
Kel-F, CTFE	198-211	199	2.3
Polymethylpentene	240	175	---
Tefzel, TFE+CE	271	200	2.0
Halon, (P)TFE	>320	260	1.5
Teflon®, FEP	252-262	204	---
Teflon®, PFA	302	260	1.5
Polycarbonate	241	121	0.7
Quartz, fused	>1665	---	0.6

reprinted with permission from Kingston, H.M.; Jassie, L.B. "Introduction to Microwave Sample Preparation" copyright 1988, American Chemical Society


CEM

Components of a Microwave System

- ◆ Magnetron
- ◆ Waveguide
- ◆ Cavity
- ◆ Circulator
- ◆ Turntable
- ◆ Wavelength Attenuator Cutoff
- ◆ Duty Cycle


CEM

Schematic of CEM Microwave Components




CEM


Heating Efficiency of Microwave Cavities



Schematic of Microwave Pattern Interaction with Pressure Vessels




Double Magnetron Design



Dr. S. Leikin
SPECTRO

CEM Microwave Components



Cycling the Magnetron

Controls the power to obtain average power level

Duty Cycle = $\frac{\text{length of time magnetron is on}}{\text{time base}}$

Time base = 60 Hz (60 times/sec)

Examples

- 500 W on a 1000 W unit (50%) with a 60 Hz time base the magnetron goes on 30 times/sec and off 30 times/sec. In 10 min at 50% power, a MW field is created 18,000 times and gives the appearance of continuous power
- 500 W on a 1000 W unit (50%) with a 1 sec time base has the magnetron on 0.5 sec and off 0.5 sec. In 10 min heating at 50%, a field is created 600 times

CEM

Calibration of a Microwave Oven

- ◆ Weigh 1 kg of ~ 23°C DI water in a plastic vessel
- ◆ Measure initial temperature of the water (T_i) to ± 0.05 °C
- ◆ Irradiate 1 kg of water for 2 minutes at three different power settings (e.g., 100%, 50% and 25%).
- ◆ Measure final temperature of water (T_f), ± 0.05 °C with stirring; use highest temperature
- ◆ Repeat twice more with new sample of room temperature water in room temperature container

- ◆ Calculate unit power according to the formula


$$\text{Power} = \frac{K \times C_p \times M \times \Delta T}{t} \quad \Delta T = (T_f - T_i)$$

for 1000 g water and 120 s, the equation becomes

$$\text{Power (watts)} = 34.85 \times \Delta T$$


CEM

Linearity of Proportional Power




CEM

Water Starting Temperature and the Measurement of Full Power



CEM

Effect of Changes in Line Voltage on Absorbed Power (measured in 1 Liter of water)



CEM

Sources of Error in Microwave Calibration

Method Dependent

- Vessel Material- heat loss, absorption of MW
- Vessel Configuration- height and diameter
- Engineering Design- cavity dimensions, exhaust fan, hot spots
- Electrical- line voltage, power supply (capacitor)

Measurement Dependent

- Temperature- accuracy and size of ΔT
- Starting Temperature of Water
- Properties of Water- dielectric of dissolved ions and heat capacity
- Irradiation Time
- Low Power- errors greater than at high wattage

CEM

Absorbed Microwave Power in Multiple Vessels

Number of Vessels	Thermister	Digital Thermometer		Mean
1	989 ± 13	963 ± 5	967 ± 2	973 ± 8
2	965 ± 12	991 ± 30	997 ± 15	984 ± 19
5	982 ± 5	995 ± 10	986 ± 17	991 ± 6
\bar{X}	979 ± 12	983 ± 17	986 ± 17	982 ± 9

- uncertainty is expressed as one standard deviation
- range of uncertainty is 0.4% to 3%
- MDS-205 at 100%, starting water temperature 23 ± 2 °C, 12/3/90

CEM

Microwave Hardware Design

Objectives of Sample Preparation


- Ability to prepare multiple samples
- Ability to monitor and control reaction conditions
- Ability to achieve rapid sample throughput
- Automation

Problems

- Uniform distribution of energy to multiple samples
- Temperature and pressure monitoring in MW
- Strong microwave-transparent vessels
- Relief mechanism must be microwave immune


CEM

Microwave Closed Vessel vs. Open Vessel (Hot Plate)



Classical Digestion Approaches Open Vessel on Hot Plate


- Easy to Use
- Low Capital Cost
- Slow
- Labor Intensive
- Easily Contaminated



Classical Digestion Approaches

Closed Vessel

- Raise reagent temperature above boiling
- Sealed environment
- Limited sample size
- Pressure build up difficult to control



Organic Samples


Matrix, Size, and Reagents

Generate high gas volume during dissolution.
Sample size helps determine method.




> 0.5 g
Open Vessel

< 0.5 g
Closed Vessel



Temperature Curve for Open Vessel and Closed Vessel Microwave Heating




12 each 120 mL PFA vessels, 19 mL HNO₃ per vessel, 648 W power



CEM

Temperature Measurement Device

Condensation on Cool Vessel Walls



© H. M. Skip Kingston



Heating of Acids in a Microwave Cavity

Thermodynamic Basis of Measurement

$$P = \frac{(K) (Cp) (m) (dT)}{t}$$

P, power absorbed by sample in watts, W (W=joule/sec)
K, conversion factor for thermochemical calories to W =4.184
Cp, heat capacity, thermal capacity, specific heat, cal/g/C
m, mass of the sample in grams, g
dT, final temperature minus initial temperature, dT=(T_f-T_i)
t, time in seconds

Estimate T_f from Estimate time to reach T_f

$$T_f = T_i + \frac{(P) (t)}{(K) (Cp) (m)} \quad t = \frac{(K) (Cp) (m) (dT)}{(P)}$$

Reprinted with permission from Kingston, H.& Jassie, L. in *Introduction to Microwave Sample Preparation: Theory and Practice*, 1988, ACS, Ch 6

Microwave Sample Preparation

How to Control The Digestion ?

Reaction Chemistry and Mechanisms
 Reagents (e.g., Acids) and Temperature
 Decomposition Rate
 Separation or Oxidation vs. Explosion
 Extent of Extraction

Distribution
 Contamination of Reaction Environment
 Loss of Analyte

New Available Mechanism
 Unique T and P Capability for Chemical Reactions
 Retention of Elements Evaporated by Microwave Only

Why these concepts are so important?

- **Temperature and pressure relationships are unique in microwave equipment**
 - They are not predictable from convection and conduction data tables
- **Boiling points of many solvents are elevated and not standard under MW conditions**
- **Non-classical relationships exist in closed vessels**
- **These unique conditions aid in specific technology and applications**

Goal:

Achieve Control of Sample Preparation/ Decomposition/Extraction/ (and post reaction sample manipulation) Using Unique Mechanisms (& apparatus)

Control of

- Reaction Chemistry
- Microwave Energy Transfer
- Microwave Reaction Mechanisms
- Equipment Configuration and Operation
- Reaction Environment

Key Unique Technical Advantages of Microwave Energy Application

- Microwave enhanced chemistry is unique
- Key technical advantages in sample preparation
 - Unique relationships
 - Chemical reaction control
 - Standardization
 - Clean chemistry and its strategic effect
 - Applications
- Microwave Technologies



© H. M. Skip Kingston

Quantum Absorption of Microwave Energy is Predictable and Controllable

Fundamental Relationship

$$P_{\text{absorbed}} = \frac{K C_p m \Delta T}{t}$$

Used to Predict Temperature

$$T_f = T_i + \frac{P_{\text{absorbed}} \cdot t}{K \cdot C_p \cdot m}$$


Deviation from Ideal Conditions

$$T_f = T_i + \frac{P_{\text{absorbed}} \cdot t}{K \cdot C_p \cdot m} - \text{Heat Loss}$$

Dependent on Equipment Configuration

© H. M. Skip Kingston


Final Microwave Heating Profile of 16 mL of Nitric and Hydrofluoric Acids (5:3 v/v).
What is happening? Why do they appear this way?



(Reprinted with permission from ACS, copyright 1988 American Chemical Society)

© H. M. Skip Kingston

Temperature and Pressure Profile of 16 mL of Pure Nitric and Hydrofluoric Acids (5:3v/v) Only
No Sample Present



(Reprinted with permission of ACS, copyright 1988 American Chemical Society)

© H. M. Skip Kingston

Example: Vapor Pressure Nitric Acid

Temperature, °C	Pressure lb/in ²	P, atm
133	180	12.3
165	380	25.9
192	630	42.9
219	995	67.7
256	1565	107
285	2245	153
313	2945	200

Reference
Journal of Research of the National Bureau of Standards, vol. 30, Feb. 1943, p 110.

© H. M. Skip Kingston

Conventional Open Vessel on a Hot Plate

At Atmospheric Pressure
HNO₃, maximum (boil point) T ~ 120 °C, w/continual heat loss

Conventional Steel Jacketed Bomb in an Oven





Bomb Reaches Thermal Equilibrium with Oven
Liquid and Gas Phases Both at Final Temperature
HNO₃ at 180 °C, Pressure ~ 40 atm

Microwave Closed Uninsulated Vessel (All PFA)

Pressure is Limited by Heat Loss in Vapor Phase
HNO₃ at 180 °C, Pressure ~ 8 atm

This Condition Results in a
"Sustained Dynamic Non-Equilibrium"

© H. M. Skip Kingston



Summary




- $PV = nRT$
 - V and R are constants
 - If either T or n are reduced, P is reduced
- Microwave systems are not in equilibrium and heat loss reduces both n (the number of gas molecules) in the gas phase and the T (temperature) of the gas phase as the gas phase collides with the energy transparent vessel walls, which are being cooled externally

© H. M. Skip Kingston

Heat Loss from the Microwave Bombs

	Temperature				
	Non-Insulated Vessel		Insulated Vessel		
	PFA Vessel (Saville TM)		High Pressure (PARR TM)		
Reagent	Inside	Outside	Inside	Outside	
H ₃ PO ₄ /HNO ₃	132°C(60 mL)	92°C (Bottom)	>200°C	47°C (Side)	
HNO ₃ /HCl	78°C(120 mL)	72°C (Bottom)	>200°C	44°C (Bot.)	

© H. M. Skip Kingston










Absorption Efficiency

Absorption Efficiency, % = $\frac{Av \text{ Power }^*, W \times 100}{\text{Input Power, W}}$

* measured in the first 90 seconds


CEM





- ### Safety in the Analytical Laboratory Using Microwave Decomposition
- Radiation**
 - Leak Detection
 - Safety Interlocks
 - Wavelength Attenuator Cutoffs
 - Remote Operation
 - Pressure**
 - Vessel Construction
 - Safety Release Devices
 - Software
 - Temperature**
 - Materials and Construction Design
 - Fiber Optic Thermometer
 - Chemistry**
 - Acid Choice
 - Matrix Composition

Sample Heating by Microwaves



Volatile Elements and MW Heating in Closed Systems

Normal Heating

- ◆ Solvated ions (Cl^- , F^-) in solution have no vapor pressure
- ◆ At elevated temperatures, vapor pressure of Cl^- and F^- metal salts is much higher than acid vapor pressures

Reduced Pressure (vacuum)

- ◆ Solution vaporization temperatures appear to decrease as volume decreases; accompanying boiling points and acid vapor pressures decrease as well
- ◆ Final solution temperature of 3 mL is $< 60^\circ\text{C}$. B.P. of volatile salt is never reached

Link, D. D., Kingston, H. M. *Anal. Chem.* 72(13), July 1, 2000, p 2908




Safety in Microwave Digestion Systems

“Laboratory Microwave Devices are Chemical Reaction Systems”

- ◆ Microwave digestion systems are general purpose systems. Reactants and reaction conditions are not specified and are often unknown in some cases.
- ◆ Microwave digestion systems are designed to meet electrical, mechanical and chemical safety standards, as well as safety factors specific to microwave heating
- ◆ Microwave digestion systems that have a means of cooling (air flow or liquid) remove heat from outer jacket and can moderate reaction rates
- ◆ Microwave digestion systems do not control pressure directly (*i.e.*, no control or check valve or back pressure regulator).



Multiwave: Temperature and Pressure Control



Mechanical and Radiation Safety

Mechanical


- Door- positive button release; interlock monitoring system
- Latches- safety interlocks; separate circuitry
- Thermal Switches- prevent magnetron overload, door interlocks compromised
- Exhaust-variable speed fan; corrosion resistant plastic hose
- Inlet/Outlet Ports- stainless steel wavelength attenuators
- Isolator-patented reflected power circulator
- Teflon Coated Stainless Steel Cavity- $< 5 \text{ mw/cm}^2$ leakage


Radiation

- Microwave Leakage Detectors- survey meters



Pressure Release Door





Pressure Can be Dangerous

Vessels

- Construction- molding, machining
- Design- threaded, pressure seals, edges
- Materials- plastic, polymers, glass, metal rotors, and frame construction

Safety Devices

- Relief disks
- Relief diaphragms, membranes
- Compression devices
- External relief valves

CEM

**MARS 5 Digestion System
Hardware & Software Safety Features**

Hardware (integral software)

- ReactiGuard™: sensor disables magnetron in case of disruptive event in the cavity
- Turntable sensor: senses stop-software override restarts rotation
- TempGuard™ (optional): IR temperature sensor to prevent vessels overheating
- SafetyLock Door: positive button release, spring-loaded metal door (burps when vessel vents violently)
- Safety Switch: shuts down magnetron in case of overheating
- Isolator: shunts reflective MW energy to dummy load, prevents magnetron from overheating

CEM

**MARS 5 Digestion System
Hardware & Software Safety Features**

Software

- Temperature: 0-300 °C (jacketed); 0-260 °C (Thermo-Optic) automatic default at 210 °C
- Pressure: sensed 200 times/minute- control to 800 psi, and monitor to 1500 psi; sensor drop > 20psig/5 sec shuts off MW power
- AutoLoad™ sensing: checks power to maintain 90% on rate

CEM

Chemical Safety Concerns at High Temperature

Mineral Acids

Perchloric.....	Dangerous hot Explosive with potassium Decomposes to Cl ₂ gas
Sulfuric.....	Dehydrating agent
Hydrofluoric..	Biological irritant/poison
Aqua Regia....	Nitrosyl chloride gas irritant

Alkaline Hydroxides

NaOH, KOH, LiOH..... Caustic, dehydrating, biological irritants

Peroxides

Hydrogen.....	Potent oxidizer
Organic Ethers..	Explosive

Organic Solvents

toxicity; explosiveness; flammability; noxiousness; volatility

CEM

Unique Hazard of Metal/Alloy Acid Digestions in a Microwave Unit

- ◆ Metals with negative reduction potentials liberate hydrogen gas; e.g., Pb²⁺, Sn, Ni, Zn, Cr, Fe²⁺, Mn²⁺, Mg, Na, Li
- ◆ Samples sealed in air
- ◆ Potentially flammable/explosive mixtures may form where metals mixtures' activation energy to ignite is very low
- ◆ Interaction of metal particles and strong magnetic field can generate sparks

Example: Titanium metal shavings, 0.1g in HCl:HF, 10:5 mL
2 ACVs with Temperature and Pressure control

Digestion program	1	2	3
power	80	80	80
pressure	100	150	200
run time	10	30	30
TAP	5	20	20
temp, °C	180	180	180

note: After 3 min, cover was blown off non-control vessel

CEM

How to Avoid Forming Potentially Dangerous Hydrogen Mixtures During Microwave Acid Digestion of Metals/Alloys

- ◆ Seal closed digestion vessels in inert gas atmosphere
 - ◆ *e.g.*, glove-box under nitrogen or argon
- ◆ Purge vessel with nitrogen or argon after addition of acid
- ◆ Purge microwave compartment with inert gas (argon or nitrogen) when open vessels are used



Compounds Unsuitable for Closed Vessel Microwave Acid Digestion

- ◆ Explosives (TNT, nitrocellulose, *etc.*)
- ◆ Propellants (hydrazine, ammonium perchlorate, *etc.*)
- ◆ Pyrophoric chemicals
- ◆ Hypergolic mixtures (nitric acid with phenol, triethylamine, or acetone)
- ◆ Animal Fats (glycerol esters undergoing nitration to nitroglycerin)
- ◆ Aviation Fuels (JP-1)
- ◆ Acetylides (compounds of acetylene)
- ◆ Glycols (ethylene glycol, propylene glycol, *etc.*)
- ◆ Perchlorates (potassium, ammonium)
- ◆ Ethers (Cellosolve, *etc.*)
- ◆ Lacquers
- ◆ Alkanes (butane, hexane, *etc.*)
- ◆ Ketones (acetone, methyl ethyl ketone, *etc.*)



Microwave Sample Preparation Techniques



Acid digestion for atomic absorption (AA), emission (ICP) and mass spectroscopy



Solvent extraction for gas and liquid chromatography



Hydrolysis of proteins and peptides for amino acid Analysis



Sample drying and moisture determination



Polymer dissolution for molecular weight determination by gel permeation chromatography (GPC)



Organic synthesis reactions



Acid digestion of reinforced composites for gravimetric determination of fiber content



Dry ashing

